

**New &
Improving**

*You're
the Artist*

SPIDENT

Variety of Indications with EsFlow

EsFlow[®]

Light Cured Nanohybrid Flowable Composite Resin

Advantages

- Highly esthetic
- Superior compressive strength
- Low polymerization shrinkage
- Flow on demand handling
- Optimal radiopacity
- Strength equivalent to natural enamel after light-cured
- Available in 2 viscosity : EsFlow[®], EsFlow[®] LV
- Available shades : EsFlow[®] : A1, A2, A3, A3.5, B1, B2, Blue
EsFlow[®] LV : A1, A2, A3



SPIDENT CO., LTD.

203 & 312, Korea Industrial Complex, 722, Gojan-Dong, Namdong-Gu, Incheon, Korea 405-821

Tel. 82-32-821-0071 Fax. 82-2-318-0034 E-mail. synergy@spident.co.kr www.spident.co.kr

CE
0120

EsFlow[®]



Light Cured Nanohybrid Flowable Composite Resin




• Indication

- Restoration in Class III~V
- Deep Class I applicable
- Base / Liner under direct restorations
- Undercut blockout
- Usable in Orthodontics
- Repair of composite
- Pit & fissure sealant

• Physical Properties




Fluidity



EsFlow (SPIDENT)

EsFlow LV (SPIDENT)

Slumping



EsFlow (SPIDENT)

Product D

Product E

Product F

• Delivery Form

Syringe Type

EsFlow[®] Syringe Refill (Pack)

Contains : 2g x 2 syringes, 8 Dispensing tips
Shade : A1, A2, A3, A3.5, B1, B2, Blue



Syringe Type

EsFlow[®] LV Syringe Refill (Pack)

Contains : 2g x 2 syringes, 8 Dispensing tips
Shade : A1, A2, A3

