

You're  
the Artist

SPIDEN<sup>T</sup>

# NOP<sup>®</sup>

Silicone Coating  
Significantly Reduces  
Pain and Tissue Trauma

## Sterile Single Use Dental Needles

- Tri-bevel point for easy tissue penetration without coring
- High tech steel tube, safe, non-toxic, pyrogenic free
- Round cutting edge protects the gum

**27G/S**  
(0.4 x 21mm)

**27G/L**  
(0.4 x 30mm)

**27G/L**  
(0.4 x 35mm)

**27G/L**  
(0.4 x 38mm)

**30G/XS**  
(0.3 x 12mm)

**30G/S**  
(0.3 x 16mm)

**30G/S**  
(0.3 x 21mm)

**30G/L**  
(0.3 x 25mm)

**31G/XS**  
(0.26 x 12mm)


CE  
0120

You're  
the Artist


SPIDENT

## Sterile Single Use Dental Needles


# NOP<sup>®</sup>



27G - 3D Cutting



30G - 3D Cutting



<b>27G/s</b> (0.4 x 21mm)	27G/Short (0.4x21mm) (blue)
<b>27G/L</b> (0.4 x 30mm)	27G/Long (0.4x30mm) (navy)
<b>27G/L</b> (0.4 x 30mm)	27G/Long (0.4x35mm) (brown)
<b>27G/L</b> (0.4 x 35mm)	27G/Long (0.4x38mm) (purple)
<b>30G/XS</b> (0.3 x 12mm)	30G/X-short (0.3x12mm) (yellow)
<b>30G/s</b> (0.3 x 16mm)	30G/Short (0.3x16mm) (pink)
<b>30G/s</b> (0.3 x 21mm)	30G/Short (0.3x21mm) (orange)
<b>30G/L</b> (0.3 x 25mm)	30G/Long (0.3x25mm) (red)
<b>31G/XS</b> (0.26 x 12mm)	31G/X-short (0.26x12mm) (green)



100 Needles



SPIDENT CO., LTD.

203 & 312, Korea Industrial Complex, 722, Gojan-Dong, Namdong-Gu, Incheon, Korea 405-821  
Tel. 82-70-5057-0827 Fax. 82-2-318-0034 E-mail. synergy@spident.co.kr / export@spident.co.kr

[www.spident.co.kr](http://www.spident.co.kr)

CE  
0120