

viOpex[®]

Temporary Root Canal Filling Material

Premixed Calcium Hydroxide Paste with Iodoform



Applications

- Apexification
- Formation of hard tissue barrier
- Check the irregular shape of root canal
- Temporary root canal filling
- Root canal filling materials for deciduous teeth
- Root resorption

Compositions

- Calcium Hydroxide
- Iodoform
- Silicon oil
- Others


Characteristics

- Excellent antibacterial effect
- Superior radiopacity
- Biocompatibility
- High flowability
- Easy to handle (syringe type)

Components

- 2.2 g Syringe, Disposable tip
- Endo-stop
- Tip cap

Radiopacity



Physical Properties

