

viOpaste®

Temporary Root Canal Filling Material

Premixed Calcium Hydroxide Paste with Barium Sulfate



Applications

- Apexification
- Formation of hard tissue barrier
- Check the irregular shape of root canal
- Temporary root canal filling materials
- Root canal filling for exudation control
- Pulp capping
- Root resorption

Compositions

- Calcium Hydroxide
- Barium Sulfate
- Hydroxyapatite
- Others


Characteristics

- Easy to remove with good water solubility
- Excellent antibacterial effect
- Superior radiopacity
- Strong bactericidal effect
- High flowability
- Excellent biocompatibility
- Easy to handle (syringe paste type)

Components

- 2.0 g Syringe, Disposable tip
- Endo-stop
- Tip cap

Radiopacity



Physical Properties

